



Gwanara, style and elegance...

LI602 Headboard

GWAST**
L: 1078
H: 712

GWANO**
L: 936
H: 543

CH700 A Bedside cabinets

2MH*
L: 500
D: 450
H: 750

2MO*
L: 400
D: 450
H: 600

4MO*
L: 400
D: 450
H: 750

5MO*
L: 500
D: 450
H: 900

5MO 2P
L: 500
D: 450
H: 900

ME710 A Chests of drawers

4MO A4*
L: 800
D: 450
H: 750

4MO A3*
L: 800
P: 450
H: 750

3MO A3*
L: 1000
D: 450
H: 750

3MO A2*
L: 1000
D: 450
H: 750

2MO A4*
L: 1200
D: 450
H: 750

2MO A3*
L: 1200
D: 450
H: 750

3MO A5
L: 800 W3
L: 1000 W4
D: 550
H: 750

CH710 A Desks

2MS A4*
L: 1400
D: 450
H: 750

1MO A2*
L: 1200
D: 450
H: 750

ME710 A Desks

4MO
L: 1740
D: 450
H: 750
Four-drawer option

CH700 A Desks

4MO
L: 1340
D: 450
H: 750

4MO A4
L: 1340
D: 450
H: 750

ME740A Wardrobe

R1
L: 800 W3
L: 1000 W4
D: 550
H: 450

W2*
L: 500
D: 550
H: 1950
Available without drawers

W3*
L: 800
D: 550
H: 1950
Available without drawers

W4*
L: 1000
D: 550
H: 1950
Available without drawers

W7*
L: 1500
D: 550
H: 1950
Available without drawers

Quality and durability



- 2 mm rounded edges and no sharp corners.
- Adjustable hinges for doors that close perfectly.
- Selection of top quality wood materials.
- 0,8 mm laminate on horizontal surfaces.
- Solid wood base.
- High quality manufacturing that ensures a high resistance of the furniture.
- Environmentally-friendly finishing.

AFNOR CERTIFICATION - 11, avenue Francis de Pressensé, 93571 SAINT-DENIS LA PLAINE cedex - France
 * NF MEDICAL "Furniture" compliant.
 Hill-Rom products are certified NF Healthcare Furniture. They are designed, manufactured and checked to offer the best possible guarantees.
 The technical reference system guarantees that the products are suited to this very specific environment:
 HYGIENE, DURABILITY, SAFETY, STRENGTH, COMFORT.
 ** Bears the NF Medical sign - Authorisation n° 178-01/01
 In accordance with the following standards : NF-EN-60601-1 (1996), NF-EN-60601-1-2 (2001), NF-EN-60601-2-38 et amdt 1 (1999), NF-S-90-312 (1984)
 Complies with essential requirements of EC directive 93/42/EEC applicable to classe 1 medical equipment
 Development and production site certified for Quality and Environmental Protection:
 ISO 9001 (2000) LNE n° 01391 - ISO 13485 (2003) G.MED n°0851 - ISO 14001 (2004) LNE n°11578-0




www.hill-rom.com

France.....	+33 (0)2 97 50 92 12	Switzerland.....	+41 (0)21/7062130	Portugal.....	+351 (0) 210 991 8910
United Kingdom.....	+44 (0)1530 41 1000	(deutschsprachig).....	021/706 21 38	Nordic Countries.....	+46 (0)8 564 353 60
Germany.....	+49 (0)211/16450-0	Austria.....	+43 (0)2243 / 28550	Export.....	Contact your local Hill-Rom distributor
Netherlands.....	+31 (0) 347 32 35 32	Ireland.....	+353 (0)1 413 6005		or contact your Area manager via website
Italy.....	+39 02-950541	Spain.....	+34 (0)936856000		or call +1 812 934 8173

Manufactured by Hill-Rom Inc. Distributed by Hill-Rom Inc.
 Hill-Rom reserves the right to change any item, specification or model without prior notice. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
 © 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Meltis, nursing home furniture



Gwanara, style and elegance...



Enhancing Outcomes for Patients and Their Caregivers...



Enhancing Outcomes for Patients and Their Caregivers...

www.hill-rom.com - Meltis Gwanara LTC - V001 - 25.05/09 - 100 ex



Gwanara, style and elegance...



① Easy to use and attractive

Improving the working conditions of caregivers is a fundamental issue.

- Gestures - Easy access to the content of the furniture - 170° hinged doors
- Posture - Easy to use and maintain - Metal drawers
- Safety - Locks can be fitted in almost all products
- Comfort - Silent-closing pads on doors and drawers for noise reduction and protection against impacts

② Maintenance and comfort

It is essential to meet the needs for hygiene and maintenance in Nursing Homes. The design and materials of the Gwanara range have been carefully chosen for greater maintainability and improved comfort for residents.

③ Mobility & well-being

Prolonging the residents' mobility is crucial to their well-being.

The numerous points of support on the furniture and the low position of the bed contribute to the enhanced independence of the residents.

Residents are able to rely on a range of handrails for support and when moving around the room. An egress bar is also available to help residents getting out of the bed.



The elegance and refinement of prestigious furnishings...

The Meltis product line benefits from all of Hill-Rom's expertise and know-how.

You can use the four separate styles to design living areas that meet your expectations.

With its classic design, noble materials and quality finishing, the Gwanara line produces an elegant and refined environment. The Gwanara line includes a collection of furniture and accessories, from head and footboards to bedside cabinets and wardrobes, that can be used to design an intimate environment for patients, visitors and caregivers.

The harmonious atmosphere produces a sense of peace of mind conducive to the administration of care, sound relationships and well-being.



①

②

③