



Meltis Cabinet CH242

At a glance :

- Integrated Minibar
- High accessibility. Door opens 110° with invisible hinges.
- Several options of configuration
- Rounded corners and edges. No sharp angles to avoid cuts and injuries.
- Great handle choice including four soft coloured handles for visual deficient or dementia resident.

Single function for double comfort

The bedside cabinet CH242 combines the comfort of an integrated minibar with the large storage capacity of the Meltis cabinets. It simplifies the work of the caregiver and it enhances the comfort and autonomy of the resident.

Meltis Cabinet CH242



The bedside cabinet CH242 with integrated minibar is an excellent solution determined to ease caregivers' daily work and to enhance comfort and autonomy of residents.

Options :


- Laminated base on 4 plastic end fitting
- 2 castors and 2 end fittings on simple base
- 4 Casters D75mm (2 front casters with-brake)
- Refrigerator - 23 litres (base)
- Refrigerator - 30 litres
- Sliding tray (with 23 litres refrigerator only) plate
- Left or right door opening
- Security lock on drawer

Technical specifications :

Bedside cabinet CH242 has an integrated minibar of 23 or 30 Litres. Box and front are made of laminated particle board with ABS edging. Top of the cabinet is made of Medium-density fiberboard (MDF) in version M2 and with laminated particle board in version M1.

There are several options to configure the furniture base including four feet in solid beech, four casters of D75mm, or a combination of two casters and two end fittings on simple base.

Weight 45 Kg
 Length 510 mm
 Depth 470 mm
 Height 898 mm
Top drawer :
 Length x depth x height 390x330x85 mm



Certification Qualité et Environnementale du site de développement et de production :
 ISO 9001 (2000) LNE n° 01391 - ISO 13485 (2003) G.MED n°0851 - ISO 14001 (1996) ECOPASS n° F-563601



Enhancing Outcomes for Patients and Their Caregivers...

France.....	+33 (0)2 97 50 92 12	Suisse/Schweiz.....	+41 (0)21/706 21 30	Nordic Countries.....	+46 (0)8 564 353 60
United Kingdom.....	+44 (0)1530 41 1000	(deutschsprachig).....021/706 21 38	Export.....	Contactez votre distributeur Hill-Rom local
Deutschland.....	+49 (0)211/16450-0	Österreich.....	+43 (0)2243 / 28550ou votre manager régional via le site web
Nederland.....	+31 (0) 347 32 35 32	Ireland.....	+353 (0)1 413 6005ou composez le +1 812 934 8173
Italia.....	+39 02-950541	Iberia.....	+34 (0)93 685 60 00		

Hill-Rom se réserve le droit de changer tout élément, spécification ou modèle sans notification préalable. La seule garantie dont répond Hill-Rom est la garantie expresse écrite accompagnant le contrat de vente ou de location de ses produits. - © 2008 Hill-Rom Services, Inc. TOUS DROITS RÉSERVÉS.

