



Koati, simplicity and performance...

LI802 Headboard

- KOANO****
 - L: 1078
 - H: 637
- KOAST****
 - L: 1088
 - H: 670

CH700 A Bedside cabinets

- 2MH***
 - L: 500
 - D: 450
 - H: 750
- 2MO***
 - L: 400
 - D: 450
 - H: 600
- 4MO***
 - L: 400
 - D: 450
 - H: 750
- 5MO***
 - L: 500
 - D: 450
 - H: 900
- 5MO 2P**
 - L: 500
 - D: 450
 - H: 900

ME710 A Chests of drawers

- 4MO A4***
 - L: 800
 - D: 450
 - H: 750
- 4MO A3***
 - L: 800
 - D: 450
 - H: 750
- 3MO A3***
 - L: 1000
 - D: 450
 - H: 750
- 3MO A2***
 - L: 1000
 - D: 450
 - H: 750
- 2MO A4***
 - L: 1200
 - D: 450
 - H: 750
- 2MO A3***
 - L: 1200
 - D: 450
 - H: 750
- 3MO A5**
 - L: 800 W3
 - L: 1000 W4
 - D: 550
 - H: 750

CH710 A Desks

- 2MS A4***
 - L: 1400
 - D: 450
 - H: 750
- 1MO A2***
 - L: 1200
 - D: 450
 - H: 750

ME710 A Desks

- 4MO**
 - L: 1740
 - D: 450
 - H: 750

Four-drawer option

CH700 A Desks

- 4MO**
 - L: 1340
 - D: 450
 - H: 750
- 4MO A4**
 - L: 1340
 - D: 450
 - H: 750

ME740A Wardrobe

- R1**
 - L: 800 W3
 - L: 1000 W4
 - D: 550
 - H: 450
- W2***
 - L: 500
 - D: 550
 - H: 1950

Available without drawers
- W3***
 - L: 800
 - D: 550
 - H: 1950

Available without drawers
- W4***
 - L: 1000
 - D: 550
 - H: 1950

Available without drawers
- W7***
 - L: 1500
 - D: 550
 - H: 1950

Available without drawers

Meltis, nursing home furniture




Koati, simplicity and performance...

Quality and durability



- 2 mm rounded edges and no sharp corners.
- Adjustable hinges for doors that close perfectly.
- Selection of top quality wood materials.
- 0,8 mm laminate on horizontal surfaces.
- Solid wood base.
- High quality manufacturing that ensures a high resistance of the furniture.
- Environmentally-friendly finishing.

AFNOR CERTIFICATION - 11, avenue Francis de Pressensé, 93571 SAINT-DENIS LA PLAINE cedex - France
 * NF MEDICAL "Furniture" compliant.
 Hill-Rom products are certified NF Healthcare Furniture. They are designed, manufactured and checked to offer the best possible guarantees.
 The technical reference system guarantees that the products are suited to this very specific environment:
 HYGIENE, DURABILITY, SAFETY, STRENGTH, COMFORT.
 ** Bears the NF Medical sign - Authorisation n° 178-01/01
 In accordance with the following standards : NF-EN-60601-1 (1996), NF-EN-60601-1-2 (2001), NF-EN-60601-2-38 et amdt 1 (1999), NF-S-90-312 (1984)
 Complies with essential requirements of EC directive 93/42/EEC applicable to classe 1 medical equipment
 Development and production site certified for Quality and Environmental Protection:
 ISO 9001 (2000) LNE n° 01391 - ISO 13485 (2003) G.MED n°0851 - ISO 14001 (2004) LNE n°11578-0



www.hill-rom.com



| | | | | | |
|---------------------|----------------------|------------------------|---------------------|-----------------------|--|
| France..... | +33 (0)2 97 50 92 12 | Switzerland..... | +41 (0)21/7062130 | Portugal..... | +351 (0) 210 991 8910 |
| United Kingdom..... | +44 (0)1530 41 1000 | (deutschsprachig)..... | 021/706 21 38 | Nordic Countries..... | +46 (0)8 564 353 60 |
| Germany..... | +49 (0)211/16450-0 | Austria..... | +43 (0)2243 / 28550 | Export..... | Contact your local Hill-Rom distributor |
| Netherlands..... | +31 (0) 347 32 35 32 | Ireland..... | +353 (0)1 413 6005 | | or contact your Area manager via website |
| Italy..... | +39 02-950541 | Spain..... | +34 (0)936856000 | | or call +1 812 934 8173 |

Manufactured by Hill-Rom Inc. Distributed by Hill-Rom Inc.
 Hill-Rom reserves the right to change any item, specification or model without prior notice. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.
 © 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.



Enhancing Outcomes for Patients and Their Caregivers...



Koati, simplicity and performance...

① Easy to use and attractive

Improving the working conditions of caregivers is a fundamental issue.

- Gestures - Easy access to the content of the furniture - 170° hinged doors
- Posture - Easy to use and maintain - Metal drawers
- Safety - Locks can be fitted in almost all products
- Comfort - Silent-closing pads on doors and drawers for noise reduction and protection against impacts

② Maintenance and comfort

It is essential to meet the needs for hygiene and maintenance in Nursing Homes. The design and materials of the Koati range have been carefully chosen for greater maintainability and improved comfort for residents.

③ Mobility & wellness

Prolonging the residents' mobility is crucial to their well-being.

The numerous points of support on the furniture and the low position of the bed contribute to the enhanced independence of the residents.

Residents are able to rely on a range of handrails for support and when moving around the room. An egress bar is also available to help residents getting out of the bed.

The minimalist and convenient Koati range creates a modern environment...

The Meltis product line benefits from all of Hill-Rom's expertise and know-how.

You can use the four separate styles to design living areas that meet your expectations.

With its clear and gentle lines and a multitude of colors, Koati produces an environment full of vitality and freshness.

The Koati line includes a collection of furniture and accessories, from head and footboards to bedside cabinets and wardrobes.

The countless possible combinations can be used to create an atmosphere that is dynamic and modern for patients, visitors and caregivers. The gentle and calm atmosphere of the living quarters creates a sense of peace and serenity that is conducive to the administration of care, sound relations and well-being.

